

FEATURE OVERVIEW

System platform		C4-FS
C-ISP platform		•
WidexLink		•
InterEar (IE)	Feature coordination	•
	System synchronisation	•
	Partner monitor	•
EcoTech III		•
Bands/Channels		
Analysis bands		15
Processing channels		15
Signal processing		
Compression	EDRC	•
	TruSound Compression	•
	TruSound Softener	•
Adaptive directionality	HD Locator	•
	Digital Pinna	•
	FreeFocus	•
Noise reduction	IE Speech Enhancer SII	•
	Speech Enhancer (SII)	-
	Noise reduction with SIS	-
	Widex Classic NR	•
	Wind noise control	•
Feedback Cancelling	IE Feedback Cancelling	•
Frequency transposition	Audibility Extender	•
Enhanced sound quality	3D TruSound	•
	100 Hz - 9.6 kHz*	•
User feedback and control		
Acoustic indicators	SmartSpeak	•
	SmartTone	•
	IE Partner Alarm	•
	Battery alarm	•
LED		•
IE volume control		RC-/M-DEX
IE program shift		•/RC-/M-DEX
Selected inputs		
IE Zen		•
IE Zen timer		•
Audio transmission		
WidexLink audio transmission		•
Telecoil		•
DAI		-
Listening programs		
Acoustic program options		5 & Zen+ with RC-/M-DEX accessory
Phone+***		•
Accessories		
RC-DEX, M-DEX and TV-DEX		•
SCOLA		-

*Receiver dependent

**Optional

***Choose between Phone+ or Zen+